

Standardization of titrants:

1. NaOH titrant is standardized against potassium hydrogen phthalate (KHPht). Quantities of KHPht corresponding to approximately 1,2, 3, 4 and 5 mmole are weighed accurately into titration vessels and allowed to dissolve in 30mL of DI water before titrating at 5mL/min against the NaOH titrant. The results are plotted, mmole KHPht on the x-axis, mL NaOH titrant consumed on the y-axis. The molarity of the titrant is calculated from the reciprocal of the gradient, as computed by regression analysis.

2. Ba(OAc)₂ titrant is standardized against a 0.1mol/L Na₂CO₃ solution prepared from anhydrous A.R. Na₂CO₃. Aliquots are prepared according to the following table:

mL 0.1mol/L Na ₂ CO ₃	mL DI water	mmole Na ₂ CO ₃
50	0	5
40	10	4
30	20	3
20	30	2
10	40	1

Solutions are titrated at 5mL/min, the results being plotted with mmole Na₂CO₃ on the x-axis, mL Ba(OAc)₂ titrant consumed on the y-axis. The molarity of the titrant is calculated from the reciprocal of the gradient, as computed by regression analysis.

Results:

Solution comprising NaHCO₃ and Na₂CO₃ with mole ratio of ~2:1, and with total carbonate concentration ~0.1mol/L as Na₂CO₃

NaHCO₃ = 5.38±0.014 g/L, n=10

Na₂CO₃ = 7.42±0.015 g/L, n=10

Total carbonate (as Na₂CO₃) = 10.81±0.012, n=10


Calculation formulae:

$$\text{NaHCO}_3 \text{ g/L} = \frac{((\text{mL NaOH} - \text{mL blank}) \times 84.01 \times M \text{ NaOH})}{\text{mL sample}}$$


$$\text{Total Carbonate g/L} = \frac{(\text{mL Ba(OAc)}_2 - \text{mL blank}) \times 105.99 \times M \text{ Ba(OAc)}_2}{\text{mL sample}}$$

$$\text{Na}_2\text{CO}_3 \text{ g/L} = \text{Total Carbonate g/L} - \frac{\text{NaHCO}_3 \text{ g/L} \times 105.99}{2 \times 84.01}$$

Thermometric Titration Plots:



Titration of NaHCO₃ component with NaOH



Sequential titration of "total carbonate" with Ba(OAc)₂

Legend:

Red = solution temperature curve

Black = second derivative curve

Example of titrant calibration: $\text{Ba}(\text{OAc})_2$

Molarity = $1/\text{gradient} =$
 $1/0.95201 = 1.0504$

Blank = y-intercept =
 0.112mL

